



Third Party Software Versions Validated with Spacelabs Flexport Products

The Spacelabs Flexport products allow presentation of non-Spacelabs medical device data as an additional parameter on your Spacelabs monitoring system. Each non-Spacelabs medical device software version uses a unique communication protocol to interface with your Spacelabs medical device. The purpose of the protocol is to present data in a specific predictable format.

The list below identifies third party software versions validated for use with the Spacelabs Flexport products. Any other third party software version may not present data appropriately. The list of validated third party software versions will be regularly updated at:

<http://www.spacelabshealthcare.com/patient-monitoring-connectivity/connectivity/flexport-interfaces>

Manufacturer	Model	Flexport	Spacelabs Verified Firmware Version
BBRAUN	BBRAUN INFUSION PUMP	90442A	400ES
CAREFUSION	Pulmonetics LTV 1200	91436A-A	05.02.06
CAREFUSION / PURITAN	PB840	90436A-15	4-070212 85-N
DRAEGER	Babylog 8000	90436A-10B	Medibus 3.00
DRAEGER	Evita2	90436A-10-32	CAP 2.00 / Dura 4.10
DRAEGER	Evita2 Dura	90436A-10E	CAP 2.00 / Dura 4.10
DRAEGER	Evita4	90436A-16	4.22 / XL5.1
DRAEGER	Ventilator	90436A-09	Irisa 9.0 / 9.5 / 9.6 Evita 9.07
EDWARDS	VIGILENCE VGS	90442A	sw 01.10
EDWARDS	VIGILANCE II	91439B-A	v00.56
EDWARDS	VIGILEO	91439B-A	sw 01.10
EDWARDS	EV1000	91443A-A	sw 01.10
GENERAL ELECTRIC	Engstrom Carestation	91436E-A	5.05
GENERAL ELECTRIC	AVANCE	91438A-A	5.0
HAMILTON	Galileo	90436A-07	3.2
HAMILTON	Raphael	90436A-07	3.2
MAQUET	Servo-I	91436B-A	2

Manufacturer	Model	Flexport	Spacelabs Verified Firmware Version
MAQUET	FLOW-I	91443B-A	FCI V0004 (protocol)
MASSIMO	RADICAL	90442A-02	Masimo SET- V4.1.0.1 Handheld - R4.2.0.1 D Station (RadNet) - R5.1.3.2
NELLCOR	N595	90442A	V 3.3.0.0
NELLCOR	N600	90442A	1.1.2.0
ORIDION	MICROCAP MICROSTREAM	90442A-02, 03	2.75
PULSION	PICCO I	90439-A	D7.0.7 (Demo mode)
PULSION	PULSION PiCCOplus	90442A	D7.0.7 (Demo mode)
PULSION	PICCO 2	91439B-A	3.1.0.8A
RADIOMETER	TCM 4/40	91439A-A	3.01
RESPIRONICS	ESPRIT	90436A-01	OTP VERSION 4.10, FLASH 31.40
SPACELABS	SIRIUS	90519-B	700900 9.62(Control Board). 700900 10.06(Front Panel)
SPACELABS	BLEASE	90519-B	700900 9.62(Control Board) 700900 9.02(Front Panel)
SPACELABS	ARKON	91443B-A	Contains a protocol check
VIASYS	Vela	90436A-18	2.02.08
VIASYS	Avea	91436C-A	3.9B

Please refer to <http://www.spacelabshealthcare.com/en/company/trademarks> for a full listing of Spacelabs Healthcare trademarks. Other brands and product names used herein are trademarks of their respective owners.