


**PAJUNK®**

***Paediatric Anaesthesia***  
*Cannula Systems for Caudal,  
Epidural, and Plexus Anaesthesia*



**Paediatric Anaesthesia**

*Fine and precise*

## The PAJUNK® Special Program for Paediatric

As manufactured of various specialty cannula, PAJUNK® offers a complete product line designed specifically for paediatric anaesthesia in

regional anaesthesia. We offer high quality cannula for caudal, spinal, epidural, and plexus anaesthesia.

### Caudal Anaesthesia

SPROTTE® 27 G



The SPROTTE® has initiated a worldwide trend to atraumatic puncture, which has proven itself in paediatric anaesthesia applications as well.

Crawford cannula 23 G, graduated, with depth stop




The Caudal cannula with its cone point geometry according to Crawford has proven to be particularly suitable in paediatric puncture.



### Epidural Anaesthesia

PAJUNK® produces two special types of cannulas for Epidural anaesthesia:


With the well proven Tuohy tip:



According to Marhofer with the well proven Tuohy tip:



The Paediatric set was developed in cooperation with Dr. Marhofer.




SPROTTE® cannula



Crawford tip



Tuohy tip




Tuohy tip


# ric Anaesthesia

## Plexus Anaesthesia


### UniPlex NanoLine cannula „Single Shot“: (UP/UPS) cannula



The UniPlex NanoLine cannulas were developed by PAJUNK® specifically for the disposable application in the peripheral conductive anaesthesia. These cannulas come in many different gages and sizes. PAJUNK® understands the growing trend of this application in paediatric anaesthesia so we offer two versions a SPROTTE® tip and a Facet tip with a shorter lengths and very small gages.




SPROTTE® tip




Facet tip

### PlexoLong NanoLine according to Meier „Catheter Technique“



The patented Plexus Catheter Set was developed in cooperation with Dr. Meier specifically for the continuous peripheral techniques for regional anaesthesia for children.



SPROTTE® SPECIAL tip



## Caudal

Product	Size	Item No.	PU
SPROTTE® cannula	27 G x 35 mm	111151-27A	25
SPROTTE® cannula	25 G x 35 mm	001151-29E	25
SPROTTE® cannula	24 G x 35 mm	001151-30E	25
SPROTTE® cannula	22 G x 50 mm	071151-30C	25

Crawford cannula	23 G x 25 mm	1149-12D025	25
Crawford cannula	22 G x 40 mm	1149-12E040	25

## Epidural

Product	Size	Catheter	X-ray contrast catheter Guide 0,2 µm flat filter Tuohy Borst adapter LOR syringe	Item No.	PU
EpiLong Set I					
Tuohy cannula	19 G x 45 mm	23 G x 50 cm	• • • •	0131153-50	10
Epidural cannula					
Tuohy cannula	18 G x 50 mm			1150-4K050	25
Tuohy cannula	19 G x 50 mm			1150-4I050	25
Paediatric set acc. Marhofer					
Tuohy NanoLine cannula	21 G x 40 mm	25 G x 50 cm	• • • •	0131156-51	10

## Plexus

Product	Size	Catheter	Item No.	PU
PlexoLong Set with SPROTTE® SPECIAL tip	21 G x 30 mm	25 G x 50 cm	541156-31B	10

## NanoLine

Product	Size	Item No.	PU
UniPlex NanoLine Facet Tip			
UP 3/25	24 G x 25 mm	001156-75	10
UP 3/40	22 G x 40 mm	001156-70	10
UP 3/50	22 G x 50 mm	001156-74	10
UniPlex NanoLine SPROTTE®			
UPS 3/40	24 G x 40 mm	001156-30G	10
UPS 3/50	25 G x 50 mm	00156-29G	10

**PAJUNK GmbH**  
 Medizintechnologie  
 Karl-Hall-Strasse 1  
 D-78187 Geisingen/Germany  
 Telefon +49 (0) 77 04/92 91-0  
 Telefax +49 (0) 77 04/92 91-6 00  
[www.pajunk.com](http://www.pajunk.com)